



EUROPEAN MEDICINES AGENCY  
SCIENCE MEDICINES HEALTH

1 13 December 2010  
2 EMA/CVMP/VICH/463/02  
3 Committee for Medicinal Products for Veterinary Use (CVMP)

4 **VICH GL34: Guideline on testing for the detection of**  
5 **mycoplasma contamination**  
6 **Draft**

Draft agreed by IWP	May 2002
Adoption by CVMP for release for consultation	15 May 2002
End of consultation (deadline for comments)	30 April 2003
Adoption by CVMP for release for further consultation <sup>1</sup>	9 December 2010
End of consultation (deadline for comments)	31 March 2011

7  
8

Comments should be provided using this [template](#). The completed comments form should be sent to [vet-guidelines@ema.europa.eu](mailto:vet-guidelines@ema.europa.eu)

<sup>1</sup> Draft VICH GL34 was published for consultation first in 2002. Following a 12 months it was agreed to suspend the consultation to wait for the testing of reference strains. The consultation is now re-opened to allow submission of any further comments.



# TESTING FOR THE DETECTION OF MYCOPLASMA CONTAMINATION

Recommended for Consultation  
at Step 4 of the VICH Process  
on 11 April 2002  
by the VICH Steering Committee

THIS GUIDELINE HAS BEEN DEVELOPED BY THE APPROPRIATE VICH EXPERT WORKING GROUP AND IS  
SUBJECT TO CONSULTATION BY THE PARTIES, IN ACCORDANCE WITH THE VICH PROCESS. AT STEP 7 OF  
THE PROCESS THE FINAL DRAFT WILL BE RECOMMENDED FOR ADOPTION TO THE REGULATORY BODIES OF  
THE EUROPEAN UNION, JAPAN AND USA.

## **1. INTRODUCTION**

### **1.1. Objective of the guideline**

This VICH (International Cooperation on Harmonization of Technical Requirements for Registration of Veterinary Medicinal Products) guideline is intended to facilitate the harmonized licensing of new products for veterinary use. It is important that biological products for veterinary use are free of contamination with Mycoplasmas to help assure consistency of production and final product safety. Mycoplasma contaminants may be introduced into cell culture and in ovo origin biological products through the master seeds, the master cell seed (stock), starting materials of animal origin, and in processing of biological materials during passage and product assembly. Therefore it is necessary to demonstrate through testing that Mycoplasmas are not present, within the limits of the test, in the final product, working seeds and cells and harvests, and starting materials such as the master seed, master cell seed, and ingredients of animal origin. This guideline establishes stages of manufacture to be tested and test procedures to detect the presence of Mycoplasma contamination. It will provide a unified standard that will facilitate the mutual acceptance of test data by the relevant regulatory authorities. Methods proven equivalent to the guideline method by scientifically accepted criteria could also be acceptable.

### **1.2. Background**

The present methods for testing for Mycoplasma contamination are described in the Japanese “Minimum requirements of biological products for animal use (1987)”, the European Pharmacopoeia (4th Edition, 2002, 2.6.7), and the United States Code of Federal Regulations, Title 9, 113.28. These requirements are all similar in that they require testing for Mycoplasma contamination using a broth and agar technique. The requirements do however differ in the specifics of these broth and agar tests as well as other alternative test methods that are required or approved for use in detecting Mycoplasma contamination.

### **1.3. Scope of guideline**

This guideline describes the manner in which tests conducted to detect the presence of Mycoplasma contamination in cell culture and in ovo origin biological products for veterinary use shall be done to assure the absence of Mycoplasma contamination. Tests on master seeds, master cell seeds (stocks), working seeds and cells, ingredients of animal origin, harvests and live final vaccine and harvests for killed products are included. Bacterial products which grow in the mycoplasma test media and products for which mycoplasma contamination risk has been addressed through a validated mycoplasma inactivation procedure will be considered outside this guideline. The absence of mycoplasma contamination in eggs used for production is controlled by appropriate testing of the flock, which is not covered by this guideline.

### **1.4. Test Methods**

The guideline describes two test methods: 1) expansion in broth culture and detection by colony formation on nutrient agar plates; and 2) expansion in cell culture and characteristic fluorescent staining of DNA (a technique capable of detecting non-cultivable strains). The use of the polymerase chain reaction (PCR) is currently under consideration by regulatory authorities for more rapid detection in broth and cell cultures, confirmation, and strain identification. Its use is encouraged in parallel testing to further develop, compare, and refine the technique for possible inclusion in future versions of this guideline.

## 2. GUIDELINE FOR TESTING FOR MYCOPLASMA CONTAMINATION

### 2.1. General test procedures for detecting Mycoplasma contamination

The culture method using broth and agar is the fundamental method of Mycoplasma detection. A solid and liquid media culture method shall be used to test harvests or final batches of vaccine, and ingredients of animal origin. Master seed, master cell seed (stock), and working seed and cell lots shall be tested using both a solid and liquid media culture method and an indicator cell culture method with DNA stain. Should either method result in a positive test for mycoplasma the sample is considered positive and is unsuitable for use.

Material	Broth & Agar Culture	DNA Stain
Master Seed & Master Cell Seed	Required	Required
Working Seed & Working Cell Seed	Required	Required
Ingredient of Animal Origin <sup>1,2</sup>	Required	
Harvest	When testing required <sup>3</sup>	
Final Product	When testing required <sup>3</sup>	

<sup>1</sup> Excluding eggs

<sup>2</sup> Unless a validated mycoplasma inactivation procedure has been applied

<sup>3</sup> The competent authorities require testing of different combinations of harvests and final product.

### 2.2. Culture test system validation

The culture method should be carried out to validate the sensitivity of a laboratory's mycoplasma detection method. A sufficient number of both solid and liquid media shall be used to insure the growth of a low level (approximately 100 CFU) of the following 5 strains of mycoplasmas.

*Acholeplasma laidlawii*

*Mycoplasma hyorhinis*

*Mycoplasma orale*

*Mycoplasma synoviae*

*Mycoplasma fermentans*

The species were selected to reflect a range (within a practical number) of antibiotic sensitivity (to detect inhibition of mycoplasma growth in the assay), fastidiousness, rapidity of growth, likelihood of being a contaminant, and pathogenicity in avian or mammalian target species. *Acholeplasma laidlawii* is a common cell culture contaminant of animal and possibly environmental origin. *Mycoplasma hyorhinis* is fastidious, is a common cell culture contaminant of animal origin, and is a mammalian pathogen. *Mycoplasma orale* is antibiotic sensitive and is a common cell culture contaminant of human origin.

*Mycoplasma synoviae* is fastidious (having a nicotiamide-adenine-dinucleotide [DPN, NAD] and cysteine requirement) and is an avian pathogen. *Mycoplasma fermentans* is a slow-growing organism and a common cell culture contaminant of human origin.

Working Reference of the strains used to validate the laboratory mycoplasma contamination culture test system should be approved for use by the competent authority, of low passage level, and identified relative to type culture isolates, (see Appendix 3.2 for information regarding the intended strategy to produce valid working references). The reference strains used to validate the culture test system will be appropriate to the products tested (see table). Validation for *M. synoviae* is required when materials of avian origin are used at any stage in development and production. Validation for *M. hyorhinis* and *A. laidlawii* is required when materials of mammalian origin are used at any stage in development and production. Validation for *M. orale* is required when an antibiotic has been used at any stage in development and production. The Working References

Preparations shall be used to validate each production lot of broth and agar. At least one working reference strain must be used as a control with each test.

**Required Reference Organisms by: product type; test method, and presence of antibiotics**

Vaccine type Antibiotic content Test Method	<i>A. laidlawii</i>	<i>M. orale</i>	<i>M. hyorhinis</i>	<i>M. synoviae</i>	<i>M. fermentans</i>
Avian <i>in ovo</i> origin vaccine Without Antibiotics Broth/Agar Method				X	X
Avian <i>in ovo</i> origin vaccine With Antibiotics Broth/Agar Method		X		X	X
Avian cell culture origin vaccine Without Antibiotics Broth/Agar Method	X			X	X
Avian cell culture origin vaccine With Antibiotics Broth/Agar Method	X	X		X	X
Mammalian cell culture origin vaccine Without Antibiotics Broth/Agar Method	X		X		X
Mammalian cell culture origin vaccine With Antibiotics Broth/Agar Method	X	X	X		X
DNA Staining Method Vaccine Without Antibiotics		X	X		
DNA Staining Method Vaccine With Antibiotics		X	X		

## 2.3. Culture Method

### 2.3.1. Incubation conditions

Incubate the broth culture medium or media in tightly stoppered containers in air. Incubate all agar plates under microaerophilic conditions (nitrogen containing 5-10% CO<sub>2</sub>). For the solid medium or media, maintain an atmosphere of adequate humidity to prevent desiccation of the agar surface.

### 2.3.2. Nutritive properties of a new batch of medium

Each new lot (batch) of medium must be tested for the nutritive properties using working references specified above in **Section 2.3**. Each testing laboratory must determine the inoculum for each of their working references that will contain a low level (approximately 100 CFU). Inoculate the solid medium with a low level (approximately 100 CFU) per 60 mm plate and per 100 ml container of broth

medium. Use at least one agar plate and broth container for each working reference. Incubate the agar and broth media and make subcultures from the broth onto agar at the specified intervals. The agar medium batch complies with the test for nutritive properties if approximately 100 CFU are achieved for all the working references specified. The broth complies if Mycoplasma growth on those agar plates subcultured from the broth is achieved for each working reference specified. Media formulations found effective are included in **Appendix 3.1** of this guideline.

### **2.3.3. Inhibitory substances**

Carry out the test for nutritive properties in the presence and absence of the material to be tested at the time of prelicense and whenever there is a change in the production method that may affect the detection of mycoplasmas. If growth of the working references occurs more than one sub-culture sooner without the test material than with the test material, or if plates directly inoculated with the test material have less than one-fifth the colonies of those directly inoculated without the test material, the test material contains inhibitory substances. These substances must be neutralized or their effect otherwise countered, e.g., through passage in substrates not containing inhibitors or dilution in a larger volume of medium, before the test for mycoplasma contamination is carried out. For the dilution technique, larger medium volumes may be used or the inoculum volume may be divided among multiple 100 ml flasks. The effectiveness of the neutralization or other process is confirmed by repeating the test for inhibitory substances after neutralization.

### **2.3.4. Test method**

**2.3.4.1** The amount of inoculum for each plate of solid medium is 0.2 ml of product to be examined. When an assay for mycoplasma concerns master and working seeds, master and working cells, and ingredients of animal origin a volume of not less than 10 ml of undiluted sample shall be tested in each liquid medium. The volume of final product to be tested in each liquid medium shall be as required by the regulatory authority issuing the marketing authorization. These are currently not less than 1 ml in Japan and the US and not less than 10 ml in the EU. Incubate the agar plates at 35°C to 38°C, microaerophilically, for 10-14 days in an atmosphere of adequate humidity to prevent desiccation of the surface. Incubate the liquid media at 35° C to 38° C in tightly stoppered containers in air for 20-21 days. At the same time incubate an uninoculated 100 ml portion of each liquid medium and agar plates as a negative control. If any significant pH change occurs upon the addition of the product to be examined (this should be determined at the time of prelicense), the liquid medium shall be restored to its original pH value by the addition of a solution of either sodium hydroxide or hydrochloric acid. On day 2 or 3 after inoculation, subculture each liquid culture by inoculating at least 1 plate of each solid medium with 0.2 ml and incubate them at 35°C to 38°C microaerophilically for 10-14 days. Repeat the procedure on the 6th or 7th day, again on the 13th or 14th day and again on the 20th or 21st day of the test. Incubate those agar plates inoculated on day 20 or 21 for 7 days. Observe the liquid medium or media every 2 or 3 days and if a color change occurs, subculture. Color change detection requires the addition of phenol red to the media.

**2.3.4.2** If the liquid medium or media shows bacterial or fungal contamination, repeat the test. If it is not possible to read at least one plate per inoculation day, the test must be repeated.

**2.3.4.3** Include in the test, positive controls prepared by inoculating a low level (approximately 100 CFU) of at least one of the working reference species onto the agar plates and into the broth medium or media. If the test is run on a routine basis, the control species should be rotated on a regular basis. This control shall be used in each test conducted with a medium that has been validated for nutritive properties using working references determined by the types of products being tested as specified in Section 2.3 of this guideline.

### **2.3.5. Judgment of the culture method**

At the end of the incubation period, examine all the inoculated solid media microscopically for the presence of mycoplasma colonies. The product is negative for Mycoplasma contamination if the growth of typical Mycoplasma colonies has not occurred on any of the inoculated solid media. If growth of typical Mycoplasma colonies has occurred on any of the solid media, the test and sample tested are considered positive for Mycoplasma contamination. The test is invalid if the positive controls do not show growth of mycoplasma on at least one subculture plate or the negative controls are positive for mycoplasma contamination. If either of the controls is invalid the test must be repeated. If suspect colonies are observed, confirmation of mycoplasma contamination may be accomplished using an appropriate and validated method.

## **2.4. Indicator cell culture method**

Cell cultures are stained with a fluorescent dye that binds to DNA. Mycoplasmas are detected by their characteristic particulate or filamentous pattern of fluorescence on the cell surface, and if contamination is heavy, in the surrounding areas. Mitochondria in the cytoplasm may be stained, but may be differentiated from mycoplasma.

### **2.4.1. Validation of the indicator cell culture method**

Using a VERO or other equivalent in efficiency indicator cell culture substrate, validate the procedure using an inoculum of a low level of appropriate working references of *M. hyorhinis* and *M. orale*. Both working references must be positive when stained with the DNA stain at the end of the test.

If for viral, etc., suspensions the interpretation of results is affected by cytopathic effects, the virus may be neutralized using a specific antiserum that has no inhibitory effects on mycoplasmas, or an alternative cell culture substrate that does not allow the growth of the virus may be used. To demonstrate the absence of inhibitory effects of serum, carry out the positive control tests in the presence of neutralizing antiserum. Antiserum lots may be qualified once rather than at use.

### **2.4.2. Test method**

**2.4.2.1** Seed the indicator cell culture at a suitable density that will yield confluence of the cells after 3 days of growth (example:  $2 \times 10^4$  to  $2 \times 10^5$  cells per ml,  $4 \times 10^3$  to  $2.5 \times 10^4$  cells/cm<sup>2</sup>) in a cell culture vessel of not less than 25 cm<sup>2</sup>. The indicator cell culture should be sub-cultured without antibiotic prior to use. Inoculate 1 ml of the sample to be examined into the cell culture vessel and incubate at 35° C to 38° C.

**2.4.2.2** After at least 3 days of incubation and the cells have grown to confluence, make a subculture onto cover slips in suitable containers or on some other surface (chambered slides) suitable for the test procedure. Seed the cells in the second subculture at a low density so that they reach only 50%



confluence after 3-5 days of incubation. Complete confluence must be avoided because it impairs visualization of mycoplasmas after staining.

**2.4.2.3** Remove medium from cover slips or chambered slides. Rinse the monolayer of indicator cells with phosphate buffered saline (PBS) and then fix with glacial acetic acid/methanol (1 to 3) or some other suitable fixing solution.

**2.4.2.4** Remove the fixing solution and discard. Wash the fixing solution with sterile water and dry slides completely if they are to be stained more than one hour later.

**2.4.2.5** Add a suitable fluorescent dye that binds to DNA such as bisbenzimidazole stain (Hoechst compound 33258, bisbenzimidazole, 5 ug/L) and allow to stain for a suitable time.

**2.4.2.6** Remove the stain and rinse the monolayer with water. Mount the cover slips if applicable and examine the slides by fluorescence (for bisbenzimidazole stain use a 330 nm/380 nm excitation filter, LP 440 nm barrier filter) at 100-400 X magnification or greater.

**2.4.2.7** Compare the microscopic appearance of the test cultures with that of the negative and working reference controls, examining for extranuclear fluorescence. Mycoplasmas produce pinpoint or filaments over the indicator cell's cytoplasm. They may also produce pinpoint and filaments in the intercellular spaces. Multiple microscopic fields as validated should be examined.

### **2.4.3 Judgment of the indicator cell culture method**

The product being examined is negative for Mycoplasma contamination if there is no evidence of pinpoint or filaments of extranuclear fluorescence. If the slides inoculated with the product contain evidence of pinpoint or extranuclear fluorescence indicative of Mycoplasma the test and sample tested are considered positive for Mycoplasma contamination. The test is invalid if the positive controls do not show the presence of the appropriate extranuclear fluorescence of the reference organisms or the negative cell controls contain extranuclear fluorescence. If either of the controls is invalid the test must be repeated.



### 3. APPENDICES

#### 3.1 Suggested broth and agar formulations

##### 9 CFR Mycoplasma Broth

Heart Infusion Broth	62.5 g
Proteose Peptone #3	25.0 g
Yeast Extract	12.5 ml
1 % Thallium Acetate	62.5 ml
1 % Tetrazolium Chloride	13.75 ml
Penicillin (100,000 units/ml)	12.5 ml
Heat inactivated Horse Serum	250 ml
H <sub>2</sub> O	2425 ml

Mix all ingredients well and adjust pH to 7.9 with 10 Normal NaOH.  
Filter sterilize through a 0.2 µ filter. Dispense into sterile test vessels.  
Add DPN/L-Cysteine solution before use, 2 ml/100 ml of broth.

##### 9 CFR Mycoplasma Agar

Heart Infusion Agar	25 g
Heart Infusion Broth	10 g
Proteose Peptone #3	10 g
1% Thallium Acetate	25 ml
H <sub>2</sub> O	995 ml
Heat Inactivated Horse Serum	126 ml
Yeast Extract	5 ml
Penicillin (100,000 units/ml)	5.2 ml
DPN/L-Cysteine	21 ml

Combine heart infusion agar, heart infusion broth, proteose peptone #3, Thallium acetate, and H<sub>2</sub>O.

Mix and bring to boil, then cool. Adjust the pH to 7.9 with 10 Normal NaOH.

Autoclave 20 min. at 121° C. Cool in water bath to 56° C.

Aseptically add: horse serum, yeast extract, Penicillin, and DPN-Cysteine.

Dispense 12 ml into each 15 X 60 mm petri dish.

##### DPN/L-Cysteine solution

Nicotiamide-adenine-dinucleotide (DPN, NAD)	5 g
Q.S. with H <sub>2</sub> O to	500 ml
L-Cysteine	5 g
Q.S. with H <sub>2</sub> O to	500 ml

Mix each chemical separately until dissolved.

Mix the two solutions and filter sterilize.

##### Japanese Liquid Medium for Mycoplasma

###### Basal Medium

50 % w/v Bovine Cardiac Muscle Extract	100 ml
Meat Peptone	10 g
Sodium Chloride	5 g
Glucose	1 g
Sodium L-glutamate	0.1 g

346	L-arginine hydrochloride	1 g
347	H <sub>2</sub> O	QS to 1000 ml
348		
349	Filter sterilize through 0.22 µ membrane filter or sterilize at 121° C for 15 min.	
350	Adjust the pH of the medium after sterilization to 7.2-7.4.	
351		
352	Additives for 77 ml of the Basal medium;	
353	Horse Serum	10 ml
354	Inactivated Porcine Serum	5 ml
355	25 % w/v Fresh Yeast Extract	5 ml
356	1 % w/v B-NAD (oxidized)	1 ml
357	1 % w/v L-cysteine HCL (1 H <sub>2</sub> O)	1 ml
358	0.2 % w/v phenol red	1 ml
359		
360	Previously filter sterilize the additives and aseptically add to the sterilized basal medium.	
361	The additives which can be sterilized by high pressure can be autoclaved. Penicillin G	
362	potassium, 500 units/ml of the medium, or Thallium acetate, 0.02 % w/v, can be added.	
363		
364	<b>Japanese Agar Medium for Mycoplasma</b>	
365	Basal Medium	78 ml
366	Agar	1 g
367		
368	Sterilize by autoclaving 121°C for 15 min.	
369		
370	Additives:	
371	Horse Serum	10 ml
372	Inactivated Porcine Serum	5 ml
373	25 % w/v fresh yeast extract	5 ml
374	1 % w/v β-NAD (oxidized)	1 ml
375	1 % w/v L-cysteine HCl (1 H <sub>2</sub> O)	1 ml
376		
377	Penicillin G potassium, 500 units per ml of medium, or thallium acetate, 0.02 % w/v can be	
378	added.	
379	Add the additives to basal/agar medium which has been liquefied by heating, and divide into	
380	sterile petri dishes, 45-55 mm. Cool and allow to solidify.	
381		
382	<b>EP Recommended media for detection of <i>M. gallisepticum</i></b>	
383	Liquid Medium:	
384	Beef Heart Infusion Broth (1)	90 ml
385	Horse Serum (unheated)	20 ml
386	Yeast Extract (250 g/L)	10 ml
387	Thallium Acetate (10 g/L solution)	1 ml
388	Phenol Red (0.6 g/L solution)	5 ml
389	Penicillin (20,000 I.U. per ml)	0.25 ml
390	Deoxyribonucleic acid (2 g/L solution)	1.2 ml
391		
392	Adjust to pH 7.8.	
393		
394	Solid Medium:	
395	Prepare as described for the liquid medium above but replace beef heart infusion broth	
396	with beef heart infusion agar containing 15 g/L of agar.	
397		
398	<b>EP Recommended Media for the detection of <i>M. synoviae</i></b>	
399	Liquid Medium:	
400	Beef Heart Infusion Broth (1)	90 ml

401	Essential Vitamins (2)	0.025 ml
402	Glucose monohydrate (500 g/L solution)	2 ml
403	Swine serum (inactivated at 56°C for 30 min.)	12 ml
404	B-Nicotinamide adenine dinucleotide (10 g/L solution)	1 ml
405	Cysteine hydrochloride (10 g/L solution)	1 ml
406	Phenol Red (0.6 g/L solution)	5 ml
407	Penicillin (20,000 I.U. per ml)	0.25 ml

408

409 Mix the solutions of B-nicotinamide adenine dinucleotide and cysteine hydrochloride and after  
410 10 min. Add the other ingredients. Adjust pH to 7.8.

411

412 Solid Medium:

413	Beef Heart Infusion Broth (1)	90 ml
414	Ionagar (3)	1.4 g

415

416 Adjust pH to 7.8, and sterilize by autoclaving, then add:

417	Essential Vitamins (2)	0.025 ml
418	Glucose monohydrate (500 g/L solution)	2 ml
419	Swine serum (unheated)	12 ml
420	β-Nicotinamide adenine dinucleotide (10g/L solution)	1 ml
421	Cysteine hydrochloride (10 g/L solutions)	1 ml
422	Phenol Red (0.6 g/L solution)	5 ml
423	Penicillin (20,000 I.U. per ml)	0.25 ml

424

## 425 EP Recommended Media for the Detection of Non-avian Mycoplasma

426 Liquid Medium:

427	Hank's Balanced Salt Solution (modified) (4)	800 ml
428	H <sub>2</sub> O	67 ml
429	Brain Heart Infusion (5)	135 ml
430	PPLO Broth	248 ml
431	Yeast Extract (170 g/L)	60 ml
432	Bacitracin	250 mg
433	Meticillin	250 mg
434	Phenol Red (5 g/L)	4.5 ml
435	Thallium Acetate (56 g/L)	3.0 ml
436	Horse Serum	165 ml
437	Swine Serum	165 ml

438

439 Adjust the pH to 7.40-7.45

440

441 Solid Medium:

442	Hank's Balanced Salt Solution (modified) (4)	200 ml
443	DEAE-dextran	200 ml
444	Ionagar (3)	15.65 g

445

446 Mix well and sterilize by autoclaving. Cool to 100° C. Add this to 1740 ml of the liquid  
447 medium described above.

448

## 449 EP Media Sub parts

450 (1) Beef Heart Infusion Broth

451	Beef Heart (for preparation of the infusion)	500 g
452	Peptone	10 g
453	Sodium Chloride	5 g
454	H <sub>2</sub> O	QS to 1000 ml

455		
456	Sterilize by autoclaving.	
457		
458	(2) Essential Vitamins	
459	Biotin	100 mg
460	Calcium pantothenate	100 mg
461	Choline chloride	100 mg
462	Folic acid	100 mg
463	<i>l</i> -Inositol	200 mg
464	Nicotinamide	100 mg
465	Pyridoxal hydrochloride	100 mg
466	Riboflavin	10 mg
467	Thiamine hydrochloride	100 mg
468	H <sub>2</sub> O	QS to 1000 ml
469		
470	(3) Ionagar	
471	A highly refined agar for use in microbiology and immunology, prepared by an ion-exchange	
472	procedure which results in a product having superior purity, clarity, and gel strength.	
473	It contains approximately:	
474	H <sub>2</sub> O	
475		12.2 %
476	Ash	1.5 %
477	Acid insoluble ash	0.2 %
478	Chlorine	0.0 %
479	Phosphate (calculated as P <sub>2</sub> O <sub>5</sub> )	0.3 %
480	Total Nitrogen	0.3 %
481	Copper	8 ppm
482	Iron	170 ppm
483	Calcium	0.28 %
484	Magnesium	0.32 %
485		
486	(4) Hank's Balanced Salt Solution (modified)	
487	Sodium chloride	6.4 g
488	Potassium chloride	0.32 g
489	Magnesium sulphate heptahydrate	0.08 g
490	Magnesium chloride hexahydrate	0.08 g
491	Calcium chloride, anhydrous	0.112 g
492	Disodium hydrogen phosphate dihydrate	0.0596 g
493	Potassium dihydrogen phosphate, anhydrous	0.048 g
494	H <sub>2</sub> O	QS to 800 ml
495		
496	(5) Brain heart infusion	
497	Calf brain infusion	200 g
498	Beef heart infusion	250 g
499	Proteose peptone	10 g
500	Glucose	2 g
501	Sodium chloride	5 g
502	Disodium hydrogen phosphate, anhydrous	2.5 g
503	H <sub>2</sub> O	QS to 1000 ml
504		
505	(6) PPLO broth	
506	Beef heart infusion	50 g
507	Peptone	10 g
508	Sodium chloride	5 g

H<sub>2</sub>O

QS to 1000 ml

### **Bisbenzimidazole stain solution for DNA Staining**

Hoechst compound 33258 (bisbenzimidazole), 5 µg per liter of buffered aqueous solution.

Note: The solution should be protected from light.

## **3.2 Reference production and use**

Standardization of testing and confidence in laboratory results between laboratories and between regions would be enhanced by the production and distribution of Master References common within or between regions. A program in progress at the time of this drafting includes the isolation of the 5 strains of Mycoplasma listed in **Section 2.3** by laboratories of the European Union and donated to the European Department of the Quality of Medicines (EDQM). EDQM will produce a sufficient quantity of these Master references, perform an intra-region validation/stability study, and distribute the material and data (including media formulations used in by study participants) to the 3 regional government laboratories (Japan, EU, and USA) of this VICH Mycoplasma working group. The regional government laboratories will then distribute these master references to those laboratories in their region wishing to validate their Mycoplasma testing systems. The references will be produced so that they contain approximately 100 CFU per specified inoculum. A group of laboratories in the three regions will standardize these references and validate the CFUs. The Master References will be distributed to government and biologics industry laboratories.

For Mycoplasma test validation each laboratory will be sent 3 vials of each reference strain depending on the types of products being tested. A different production lot of the media or medias shall be used for each vial of a reference strain. After the laboratory completes the testing they shall report their results to the regional government laboratory supplying the reference vials. At the time of the validation testing each laboratory shall produce and validate (including stability) working references from the Master references. These references will be produced within the passage level specified for each Master reference. These working references shall be used to test subsequent batches of Mycoplasma media used in the laboratories testing for Mycoplasma contamination in veterinary biologics.

For DNA staining validation, the following strains may also prove useful:

*M. hyorhinis* -- ATCC 29052

*M. orale* -- ATCC 23714

## **3.3. Glossary**

### **Batch (lot, serial) of starting material of animal origin**

The total quantity of homogenous material (e.g., cells, serum) identified by a unique serial number.

### **Cell-seed system**

A system whereby successive final lots (batches) of a product are manufactured by culture in cells derived from the same master cell seed. A number of containers from the master cell seed are used to prepare a working cell seed.

### **Cell lines**

Cultures of cells >10 passages or subcultures from the tissue of origin and having a high capacity for multiplication *in-vitro*.

### **Final product, batch, lot, or serial**

A collection of closed, final containers or other final dosage units that are expected to be homogeneous and equivalent with respect to risk of contamination during filling or preparation of the final product. The dosage units are filled, or otherwise prepared, from the same final bulk product, freeze-dried together (if applicable) and closed in one continuous working session. They bear a distinctive number or code identifying the final lot (batch, serial). Where a final bulk product is filled and/or freeze-dried on several separate sessions, there results a related set of final lots (batches, serials) that are usually identified by the use of a common part in the distinctive number or code; these related final lots (batches, serials) are sometimes referred to as sub-batches, subserials, sub-lots or filling lots. For the purposes of mycoplasma testing, a single sub-batch may be considered representative of the batch.

#### **Harvests**

Material derived on one or more occasions from a single production culture inoculated with the same working seed lot (single harvest) or pooled material containing a single strain or type of micro-organism or antigen and derived from a number of eggs, cell culture containers, etc. that are processed at the same time (monovalent pooled harvest).

#### **Master cell seed (stock)**

A collection of aliquots of cells (primary or cell line) of a single passage level for use in the preparation of the product, distributed into containers in a single operation, processed together and stored in such a manner as to ensure uniformity and stability and to prevent contamination. Master cell seed is usually stored at temperatures of - 70°C or lower.

#### **Master References**

The reference strains of Mycoplasma provided by regulatory authorities to satisfy the reference requirements specified in this document.

587 **Master seed**

588 A collection of closed containers of a culture of micro-organisms of a single passage level  
589 used for the production of all batches of a designated veterinary biological product,  
590 distributed from a single bulk into containers and processed together in a single operation  
591 in such a manner as to ensure uniformity and stability and to prevent contamination.  
592

593 **Microaerophilic condition**

594 A nitrogen atmosphere containing 5-10% carbon dioxide and sufficient humidity to prevent  
595 drying of the agar plates.  
596

597 **Passage**

598 One transfer of cells or microorganisms followed by the normally used incubation period  
599 for the cell or microorganism concerned.  
600

601 **Primary cell cultures**

602 Primary cell cultures are cultures of cells essentially unchanged from those in the animal  
603 tissues from which they have been prepared and being no more than 10 *in-vitro* passages  
604 to the test level from the initial preparation from the animal tissue. The first *in-vitro*  
605 cultivation is regarded as the first passage of the cells  
606

607 **Seed-lot system:**

608 A system in which successive batches of a product are derived from the same master  
609 seed virus. For routine production, a working seed virus may be prepared from the master  
610 seed virus.  
611

612 **Working cell seed (stock)**

613 A collection of aliquots of cells derived from the master cell seed and at the passage level  
614 used in the preparation of production cell cultures. The working cell seed is distributed into  
615 containers, processed and stored as described for master cell seed. The term includes  
616 production cell seed.  
617

618 **Working References**

619 A passage of the Reference strains of Mycoplasma produced in the testing laboratory for  
620 use as controls to satisfy the reference requirements specified in this document.  
621

622 **Working seed**

623 A collection of aliquots of a microorganism derived from the master seed virus and at the  
624 passage level used in the preparation of product. Working seed virus is distributed into  
625 containers and stored as described for master seed virus. The term includes production  
626 seed.