

Considerations Regarding Assessment of A Modified Monoclonal Antibody (mAb) Product Related to A Prototype mAb Product in Addressing Emerging SARS-COV-2 Variants – A CMC Perspective

Qing (Joanna) Zhou
Supervisory Biologist
Office of Biotechnology Products
Center for Drug Evaluation and Research
Food and Drug Administration

12/15/2022

Disclaimer




The views expressed in this presentation are those of the speaker and should not be used in place of regulations, published FDA guidances or discussion with the Agency.

Prototype vs Modified COVID-19 mAbs

- **Prototype mAb:** Existing COVID-19 mAbs product demonstrated to be safe and effective against previous and/or present SARS-COV 2 variants based on a clinical endpoint in randomized, controlled trial(s).
- **Modified mAb:** A new COVID-19 mAb product that is manufactured using the same platform process as the prototype COVID-19 mAb product with the main structural and functional differences limited to variable region sequences and the epitopes to which it binds.

General Principle: Leveraging Existing Knowledge





CMC Tools to Potentially Boost the Confidence in the Relevance Level



Platform Technology

- mAb construction and engineering
- mAb expression system
- mAb manufacture



Product understanding linking critical quality attributes (CQAs) to clinical performance



Product characterization

- Structure
- Function

Several CMC Points to Consider

- Potential to Reduce Uncertainty



- Platform expression system: same vectors, same methodology, same engineering/mutation (if present), differences in the construct for variable region sequences designed for the modified mAb
- Platform manufacturing process: Most of the manufacturing process and product controls for the modified COVID-19 mAb should be the “same” (e.g., unit operations, etc.) as that of the prototype COVID-19 mAb. Updated CMC data (e.g., formulation change, differences in scale and process controls, etc.) that are needed to support the manufacture of the modified COVID-19 mAb, will need to be evaluated as a scientific matter.

Several CMC Points to Consider

- Potential to Reduce Uncertainty



- Product characteristics and profile: Thorough product characterization data should be provided to demonstrate the relevance of the prototype mAb product(s) and the modified mAb product. For example,
 - Primary amino acid sequence
 - Secondary, tertiary and higher order structure
 - Post-translational modifications (deamidation, oxidation, glycosylation, etc.)
 - Biologic activities that contribute to MOA (binding affinity and kinetics, neutralization, etc.)
 - Fc effector functions including ADCC, CDC, ADCP, and other assessments
- Any specific Fc-engineering for the purpose of modifying Fc-effector functions or antibody half-life should be the same in the prototype mAb product and the modified mAb product.

Acknowledgement



- Marjorie Shapiro
- Steven Kozlowski
- DAV and OBP Staff



U.S. FOOD & DRUG
ADMINISTRATION