



Ncardia
Stem cell experts

MEA assays using human iPSC-derived cardiomyocytes; challenges and opportunities

Event


EMA Workshop 2017

Presenter

Tessa de Korte, MSc

Date

2017, October 5



Ncardia at a glance



✓	Foundation	Ncardia was established in September 2017 following a merger of Pluriomics and Axiogenesis
✓	Based in	Leiden (NL), Gosselies (BE), Cologne (DE)
✓	Team	~65 people
✓	Technology	Human iPSC-derived technology
✓	Business	<u>Products</u> – Human iPSC-derived cardiac and neural products <u>Services</u> – <i>In vitro</i> assay services for safety screening and drug efficacy screening
✓	Active in geography	Worldwide, direct presence in EU and US, distributors in Japan
✓	Clients	Pharma, Biotech, Academics labs



Chapter 1

hiPSC-derived cardiomyocytes

- Current use, challenges and opportunities

Chapter 2


CiPA

- An opportunity for hiPSC-derived cardiomyocytes: regulatory guidance for novel cardiac safety assessment

Chapter 3

Future prospects


- Moving forward with additional hiPSC-derived cell types for safety assessment




1

hiPSC-derived cardiomyocytes

**Current use, challenges
and opportunities**



Human cellular models for drug discovery and development




Human iPSC-derived cardiomyocytes

Fully functional human iPSC-derived cardiomyocytes


- obtained through *in vitro* differentiation of transgenic hiPSCs and puromycin selection technology
- obtained without any genetically modification or purification/selection procedures
- ✓ Showing the essential human physiological conditions and pharmacological responses (3R's benefit)




Pluricyte® Cardiomyocyte Kit
Cor.4U® Cardiomyocyte Kit
CardioPlate™



Alpha actin (green)
Myosin heavy chain 7 (red)





Movie of
beating monolayer



Acute pharmacology: gold standard compound E-4031 induces arrhythmia as expected

aptuit

Current clamp mode - hERG blocker E-4031



Challenges of human iPSC-derived cardiomyocytes for cardiac safety assessment



- Passive acceptance of implementation
- Large data sets preferred to better understand the benefits
- Not a multi-organ system
- Resembles the cellular electrophysiology more than human heart
- Maturation stage


- ✓ Modulating form/stiffness, investigating various chemical compounds → efficient way to provide "adult" **maturation**
- ✓ Find out the relevant parameters that could accurately **predict** different forms of cardiotoxicity (e.g. FPD ↔ QT using MEA analysis)
- ✓ Early stage drug screening
- ✓ Personalized safety testing (investigate potential disease-causing genetic variants)



Multielectrode array (MEA) assays for predicting drug-induced ion-channel blockade

Robust de- and repolarization peaks, enabling easy analysis of electrophysiological parameters:

- Field potential duration (FPD)
- Beat period
- Depolarization peak amplitude






Data generated using Axion BioSystems' Maestro MEA system

Multielectrode array (MEA) assays for predicting drug-induced ion-channel blockade


Robust de- and repolarization peaks, enabling easy analysis of electrophysiological parameters:

- Field potential duration (FPD) >> hERG channel blockade, TdP-like arrhythmia
- Beat period >> chronotropic effects
- Depolarization peak amplitude >> Na⁺ channel blockade



Data generated using Axion BioSystems' Maestro MEA system

Reproducible Field Potential Waveforms from 8 Individual Electrodes in 1 Well




Data generated using Axion BioSystems' Maestro MEA system


2

CiPA

**An opportunity for
hiPSC-derived
cardiomyocytes:
regulatory guidance
for novel cardiac
safety assessment**



Non-clinical testing strategy, ICH S7B



ICH S7B, step 4, May 2005

ICH S7B Approach

ICH S7B Approach “Successful” → No drugs removed from market due to TdP since 2003

HOWEVER:


- hERG not very predictive of QTc effects
- Findings of hERG/iKr block in early screening of drug candidates often discourages compound progression, “throwing the baby out with the bathwater”
- Some drugs might not be available today if hERG assay results “enforced” (**pentobarbital, verapamil, ranolazine**)

FINALLY:


- Translation of current *in vitro* assays to human QT findings imperfect, QT only accounts for a minor form of cardiotoxicity

CiPA: Comprehensive *in vitro* Proarrhythmia Assay, the four components



Drug effects on multiple human ion-channels



In silico reconstruction cellular human ventricular electrophysiology



In vitro Effects human iPSC-derived cardiomyocytes



Clinical evaluation unanticipated electrophysiology




Characterize/Classify Effects

Check for Missed or Unanticipated Effects




CiPA

FDA drives: *in vitro* cell based tests for the prediction of cardiac safety




ICH S7B, step 4, May 2005

Proposed modified scheme,
presented by Dr. Gary Gintant
SPS Meeting 2017




3 independent testing sites confirmed the predictivity of Ncardia's hiPSC-derived Cardiomyocytes as a human cellular model

Electrophysiological analysis of Cor.4U® on the MEA System



Analysis of different parameters by testing site

Pharmacological analysis of Cor.4U® on the MEA System




Comparison of different testing sites per compound

3

Future prospects

**Moving forward with
additional hiPSC-
derived cell types for
safety assessment**



HESI Contractility consortium using human iPSC-derived cardiomyocytes

- **GOAL**


Develop a pilot program that utilizes cellular systems for the prediction of drug-induced depression of cardiac contractility that translates to clinical cardiotoxicity.



*Human iPSC-derived CMs are a **rapidly evolving technology** that is already being evaluated for cardiac electrophysiological toxicity in the **CiPA initiative**. Their use for evaluation of **cardiac contractility** is somewhat more preliminary in development, however this methodology is also moving forward in a very promising manner, and an effort to advance that evaluation is beginning in the **Division of Applied Regulatory Science in the Office of Clinical Pharmacology at FDA**.*


Neuronal cell products follow the Cardiac products and are now considered for regulatory purposes

Peri.4U™ - Peripheral neurons




To be used for seizure liability testing
(initiated by the FDA)

Astro.4U™ - Astrocytes



BOTOX safety testing

CNS.4U™ - CNS neurons



Contact us!

www.ncardia.com



We are shaping the
future of drug discovery
& development

Tessa de Korte
Application
Specialist

Ralf Kettenhofen
Head of Product
management

Dietmar Hess
Product Manager
Neural Cells

Greg Luerman
Technical Director,
North America

Chris Fleming
Senior Application Specialist,
North America



Thank you!

Tessa de Korte, MSc
Application Specialist

tessa.dekorte@ncardia.com

