

ATTACHING NEEDLE, REMOVING AIR FROM CARTRIDGE, AND PREPARING PEN FOR INJECTION (STEPS 5-7)

- **What if I leave my pen with the needle attached (Step 5), and then come back to finish the next step much later?**

This can cause your needle to block, you should continue from **Step 5** to **Step 6** right away.

- **What if I do not attach the needle at Step 5?**

If the needle is attached **before Step 5**, some of the medicine may be lost during mixing.

DO NOT attach the needle before **Step 5**.

DO NOT turn the cartridge to **[3]** (Step 6b) before attaching the needle (Step 5).

The needle must be attached to allow air to escape while turning the cartridge from **[2]** to **[3]** in the number window.

- **What if I do not attach the needle as instructed?**

If you do not push the needle on completely, the pen may jam or leak.

If the pen is jammed or leaking, throw it away and use another pen.

- **What if I do not hear the “click” when the [2] or when the [3] is moved into the Number Window?**

If you do not hear a “click” when the **[2]** or when the **[3]** is moved into the number window, you may not have the number fully centered in the window. **Twist the clear cartridge** slightly in the direction of the arrows (clockwise) to complete the “click” and center the number in the window.

If you are unable to turn to position **[3]**, throw it away and use another pen.

REMOVING BOTH NEEDLE CAPS AND INJECTING YOUR MEDICINE (STEPS 7-8)

- **After I turn the pen to Number 3 (Step 6), there are still some small air bubbles remaining. Can I still use the pen?**

Seeing small air bubbles remaining is normal and you can still use the pen.

- **After I give my medicine, there is some liquid still seen in the clear cartridge.**

This is normal. If you have heard and felt the injection button “click” and slowly counted to **5** before pulling the needle out of your skin, you should have received the full dose of your medicine.