

<u>MA (EU) number</u>	<u>(Invented) name</u>	<u>Strength</u>	<u>Pharmaceutical form</u>	<u>Route of administration</u>	<u>Immediate Packaging</u>	<u>Content (concentration)</u>	<u>Pack size</u>
EU/1/12/807/001	Tresiba	100 U/ml	Solution for injection	Subcutaneous use	Pre-filled pen (glass) (FlexTouch)	3 ml	1 pre-filled pen
EU/1/12/807/002	Tresiba	100 U/ml	Solution for injection	Subcutaneous use	Pre-filled pen (glass) (FlexTouch)	3 ml	1 pre-filled pen + 7 NovoFine needles
EU/1/12/807/003	Tresiba	100 U/ml	Solution for injection	Subcutaneous use	Pre-filled pen (glass) (FlexTouch)	3 ml	1 pre-filled pen + 7 NovoTwist needles
EU/1/12/807/004	Tresiba	100 U/ml	Solution for injection	Subcutaneous use	Pre-filled pen (glass) (FlexTouch)	3 ml	5 pre-filled pens
EU/1/12/807/005	Tresiba	100 U/ml	Solution for injection	Subcutaneous use	Pre-filled pen (glass) (FlexTouch)	3 ml	10 (2 x 5) pre-filled pens (multipack)
EU/1/12/807/006	Tresiba	200 U/ml	Solution for injection	Subcutaneous use	Pre-filled pen (glass) (FlexTouch)	3 ml	1 pre-filled pen + 7 NovoFine needles
EU/1/12/807/007	Tresiba	100 U/ml	Solution for injection	Subcutaneous use	Cartridge (glass) (Penfill)	3 ml	5 cartridges
EU/1/12/807/008	Tresiba	100 U/ml	Solution for injection	Subcutaneous use	Cartridge (glass) (Penfill)	3 ml	10 cartridges
EU/1/12/807/009	Tresiba	200 U/ml	Solution for injection	Subcutaneous use	Pre-filled pen (glass) (FlexTouch)	3 ml	1 pre-filled pen
EU/1/12/807/010	Tresiba	200 U/ml	Solution for injection	Subcutaneous use	Pre-filled pen (glass) (FlexTouch)	3 ml	1 pre-filled pen + 7 NovoTwist needles
EU/1/12/807/012	Tresiba	200 U/ml	Solution for injection	Subcutaneous use	Pre-filled pen (glass) (FlexTouch)	3 ml	2 pre-filled pens
EU/1/12/807/013	Tresiba	200 U/ml	Solution for injection	Subcutaneous use	Pre-filled pen (glass) (FlexTouch)	3 ml	3 pre-filled pens
EU/1/12/807/015	Tresiba	200 U/ml	Solution for injection	Subcutaneous use	Pre-filled pen (glass) (FlexTouch)	3 ml	6 (2 x 3) pre-filled pens (multipack)
EU/1/12/807/016	Tresiba	200 U/ml	Solution for injection	Subcutaneous use	Pre-filled pen (glass) (FlexTouch)	3 ml	5 pre-filled pens

