

<u>MA (EU) number</u>	<u>(Invented) name</u>	<u>Strength</u>	<u>Pharmaceutical Form</u>	<u>Route of Administration</u>	<u>Immediate Packaging</u>	<u>Content (concentration)</u>	<u>Pack size</u>
EU/1/20/1513/001	YUFLYMA	40 mg	Solution for injection	Subcutaneous use	pre-filled syringe (glass)	0.4 ml (40 mg/0.4ml)	1 pre-filled syringe + 2 alcohol pads
EU/1/20/1513/002	YUFLYMA	40 mg	Solution for injection	Subcutaneous use	pre-filled syringe (glass)	0.4 ml (40 mg/0.4 ml)	2 pre-filled syringes + 2 alcohol pads
EU/1/20/1513/003	YUFLYMA	40 mg	Solution for injection	Subcutaneous use	pre-filled syringe (glass)	0.4 ml (40 mg/0.4 ml)	4 pre-filled syringes + 4 alcohol pads
EU/1/20/1513/004	YUFLYMA	40 mg	Solution for injection	Subcutaneous use	pre-filled syringe (glass)	0.4 ml (40 mg/0.4 ml)	6 pre-filled syringes + 6 alcohol pads
EU/1/20/1513/005	YUFLYMA	40 mg	Solution for injection	Subcutaneous use	pre-filled syringe (glass)	0.4 ml (40 mg/0.4 ml)	1 pre-filled syringe with needle guard + 2 alcohol pads
EU/1/20/1513/006	YUFLYMA	40 mg	Solution for injection	Subcutaneous use	pre-filled syringe (glass)	0.4 ml (40 mg/0.4 ml)	2 pre-filled syringes with needle guard + 2 alcohol pads
EU/1/20/1513/007	YUFLYMA	40 mg	Solution for injection	Subcutaneous use	pre-filled syringe (glass)	0.4 ml (40 mg/0.4 ml)	4 pre-filled syringes with needle guard + 4 alcohol pads
EU/1/20/1513/008	YUFLYMA	40 mg	Solution for injection	Subcutaneous use	pre-filled syringe (glass)	0.4 ml (40 mg/0.4 ml)	6 pre-filled syringes with needle guard + 6 alcohol pads
EU/1/20/1513/009	YUFLYMA	40 mg	Solution for injection	Subcutaneous use	pre-filled pen (glass)	0.4 ml (40 mg/0.4 ml)	1 pre-filled pen + 2 alcohol pads
EU/1/20/1513/010	YUFLYMA	40 mg	Solution for injection	Subcutaneous use	pre-filled pen (glass)	0.4 ml (40 mg/0.4 ml)	2 pre-filled pens + 2 alcohol pads
EU/1/20/1513/011	YUFLYMA	40 mg	Solution for injection	Subcutaneous use	pre-filled pen (glass)	0.4 ml (40 mg/0.4 ml)	4 pre-filled pens + 4 alcohol pads
EU/1/20/1513/012	YUFLYMA	40 mg	Solution for injection	Subcutaneous use	pre-filled pen (glass)	0.4 ml (40 mg/0.4 ml)	6 pre-filled pens + 6 alcohol pads